

CryoFMRTM

- ▶ Revolutionary variable temperature turn-key FMR spectrometer
- ▶ Extracts the temperature dependent M_S , damping, inhomogeneous broadening and gyromagnetic ratio
- ▶ Available for the Montana Instruments Cryostation the Quantum Design PPMS
- ▶ Frequency 2 – 17 GHz
- ▶ Temperatures from 10 K to 350 K
- ▶ Easy to use LabView interface

MONTANA INSTRUMENTS
COLD SCIENCE MADE SIMPLE



Quantum Design

Nanosc Instruments

info@nanosc.se www.nanosc.se

